# MEDISANA®

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# IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.

Should you give this device to another person, it is vital that you also pass on these instructions for use.

## 1 Safety Information!

- **Only use this unit according to the intended purpose and as directed in the instruction manual. The device is not intended for commercial use.**
- Before connecting the device to your power supply, please ensure that the supply voltage stated on the rating plate is compatible with your mains supply.
- Do not use any auxiliary parts that have not been recommended by the maker.
- The warranty is null and void if it is not used for the intended purpose.
- This device is not designed to be used by persons (including children) with limited physical, sensory or mental abilities, or by persons with insufficient experience and/or knowledge, unless under observation by a person responsible for their safety, or unless they have been instructed in the use of the device.
- Children must be supervised to ensure that they do not play with the device.
- When using the appliance on the skin, it might become very hot. To prevent burns to the skin, check the temperature of the attachment regularly. This applies particularly to persons who are insensitive to heat.
- Do not attempt to repair the unit yourself in the event of any defects. Not only does this render the warranty null and void but you can also run the risk of serious hazards (fire, electric shock, injury). Only have repairs carried out by authorised service points.

### ***Please note the following in particular concerning electrical equipment:***

- Only use the unit in enclosed spaces!
- Do not reach for the unit if it has fallen into water. Disconnect the mains adapter from the mains socket immediately.
- Do not use any electrical equipment while in the bath or shower. Do not store or keep any electrical equipment in places where it could fall into the bathtub or washbasin.
- Keep the mains cable away from hot surfaces and naked flames.
- Never use the unit if the power cable has been damaged, the mains adapter is defective, if the unit does not work perfectly, if it has fallen into water or has been damaged.
- Never carry, pull or rotate the unit by the mains cable and never pinch the cable between other objects.
- Uninterrupted treatment should not last longer than 15 minutes. Then switch off the unit for at least half an hour and allow it to cool down.
- Only pull the plug out of the socket when you have dry hands, when the unit has been switched off and is not being used any longer, and also do this when you have to replace or clean the heads.

### ***Important recommendations concerning your health:***

- Consult your doctor first before using this unit if you have diabetes or other related illnesses.
- Pregnant women should note the important precautionary measures and the amount of stress they can take, consult your doctor first if necessary.
- Do not treat any parts of the body that have swellings, burns, inflammation, peeling skin, wounds or other sensitive points.
- The treatment should be pleasant. Stop at once and consult your doctor if you feel pain or the treatment seems unpleasant.
- For hygiene reasons, the attachments should only be used by one person.

## Thank you for your confidence in us and congratulations!

You have acquired a quality product with your purchase of the manicure/pedicure unit from **MEDISANA**. To ensure that your **MEDISANA MP 820** provides successful results for a long time, we recommend that you read carefully through the following notes on its use and care.

## 2 Items supplied and packaging

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use it and contact your dealer or your service centre.

The following parts are included:

- **MP 820** -device
- mains adapter
- 4 accessory heads
- storage bag
- instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

### WARNING

Please ensure that polythene packing is kept away from the reach of children! Risk of suffocation!

## 3 Special points concerning the MEDISANA MP 820

- The **MEDISANA MP 820** manicure and pedicure device is simple to operate and ideal for the extensive care of nails, cuticles and slight calluses on the hands and feet. It is particularly suitable for treating ingrown and very tough nails and also for treating calluses, rough skin and corns.
- The sapphire heads have an extremely long life and virtually never wear out.
- The unit is fitted with a safe quick-clamping chuck with automatic locking.
- Do not soften the skin on hands and feet prior to treatment.
- Due to the minimal risk of injury, the **MP 820** is also especially suitable for use by diabetics or older people.
- **Why is the MP 820 especially suitable for use with diabetics?**

Because the speed is continuously adjustable and the device stops if too much pressure is used, the unit can also be used safely by persons who have reduced sensitivity to pain or are unable to react quickly to pain. Also, the sapphire accessory heads only remove dead skin.

## 4 Operating the unit

- The appliance is fitted with a NiMH rechargeable battery. You must fully charge the device before using it for the first time. Plug the mains adapter **5** into the socket at the end of the appliance and then plug the mains adapter into the power socket. The red battery charging LED **3** flashing slowly (about once every two seconds) indicates that charging is in progress. When charging is completed, the flashing stops. The LED goes out.
- Push the required attachment into the opening of the quick-clamping chuck **1** until it meets with resistance.
- The ON/OFF/speed controller **+/- 2** activates the appliance. When you slide the controller **+/- 2** about 1mm towards "+", the attachment starts turning at the slowest speed and the two LEDs **6** to the side of the attachment illuminate the area where it is being used.
- Use the switch **4** to select clockwise or anticlockwise rotation and the controller **+/- 2** to regulate the speed. Start slowly for each application and increase the continuously adjustable rotating speed to suit your individual requirements.
- This produces a small amount of dust. The sapphire heads only remove dead skin particles with careful grinding. Stop the treatment at once if you find that the removal of hard skin is unpleasant.
- If the battery is low, the performance of the device will decrease. If the battery is empty, the U-shaped LED on the head of the device is blinking. If the battery is low, mains operation is also possible.

## 5 Attachments *(supplied)*

### 1. *Sapphire cone, coarse*

For slightly callused and rough skin on feet and hands but also for abrading very thick, hard nails. It re-shapes crooked nails and relieves painful pressure.

### 2. *Sapphire disc, coarse*

Shortens and shapes nails right into the corners. The result is perfectly cared for nails which no longer split, break or peel.

### 3. *Flame file*

or careful release of ingrown nails. It is also possible to remove fungal spores from under the front edge of the nail. **Use the flame file with great care!**

### 4. *Felt cone*

The attachment for final finishing of all nail care. Finger and toe nails are buffed to a natural shine. Superbly smooth nails and nail edges round off a perfect treatment.

- To change the attachments, first pull off the part that is not required then push the desired attachment in carefully until it meets resistance.

## 6 Replacement set, article no. 85450 *(complete set of 10, optionally available)*

### *Sapphire cone, long/coarse*

For large-surface abrasion of thick callus formation and rough skin

### *Saphirscheibe, fein*

For shortening the nails by filing

### *Sapphire file, round*

For grinding out ingrown nails

### *Needle file*

For getting rid of difficult corners and edges

### *Cylindrical file*

For treating thickened, tough nails.

### *Drainage grinder*

For treating extensive, thick callus on heels and balls

## 7 Cleaning and maintenance

- Always pull the mains adapter out of the socket in the case of an electrical unit.
- Never use strong cleaning agents or stiff brushes.
- Clean the unit with a soft cloth which has been lightly moistened with a mild soap. Never dip the unit into water for cleaning and ensure that no water gets into the unit at any time. Only use the unit again, when it is fully dry.
- For hygiene reasons, all attachments may be cleaned with a brush (e.g. old toothbrush). It is also advisable to disinfect the grinding surfaces with commercially available spray disinfectant or 90% alcohol to prevent the transmission of infections. The pins of the attachments should be kept free from grease. First allow the heads to dry off completely before you use them again on the unit.
- Store the device in a dry, cool place.

## 8 Disposal

This product must not be disposed of together with domestic waste.

All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

## 9 Technical specifications

Name and model	: <b>MEDISANA</b> Manicure and pedicure device <b>MP 820</b>
Current supply:	NiMH rechargeable battery 3,6 V = ; 450 mAh (fitted in the device) Mains adapter input: 100-240 V~ ; 50/60 Hz ; 130 mA output: 5,0 V = ; 500 mA
Revolutions	: appr. 1 200 – 11 000 rpm idling
Dimensions (Hand-held unit)	: appr. 170 x 40 x 35 mm (L x W x H)
Weight (Hand-held unit)	: appr. 120 g without mains adapter
Article no.	: 85150 EAN-Number : 40 15588 85150 6
Article no. Replacement set	: 85450 EAN-Number : 40 15588 85450 7

**In accordance with our policy of continual product improvement, we reserve the right to make technical and optical changes without notice.**

## 10 Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
  - b. All damage which is due to repairs or tampering by the customer or un-authorised third parties.
  - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
  - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

**MEDISANA AG**  
Jagenbergstraße 19  
41468 NEUSS  
GERMANY

E-Mail: [info@medisana.de](mailto:info@medisana.de)  
Internet: [www.medisana.de](http://www.medisana.de)

The service centre address is shown on the attached leaflet.

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